

PPAR  $\alpha$  (Phospho-Ser21) Antibody

Catalog No: #12432

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	PPAR $\alpha$ (Phospho-Ser21) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	PPAR $\alpha$ (Phospho-Ser21) Antibody detects endogenous levels of PPAR $\alpha$ only when phosphorylated at Ser21
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human PPAR $\alpha$ (Phospho-Ser21)
Target Name	PPAR $\alpha$
Modification	Phospho
Other Names	PPARA, PPAR Alpha, PPARalpha, PPAR-alpha, NR1C1, HPPAR, PPAR
Accession No.	Swiss-Prot#: Q07869NCBI Gene ID: 5465
Uniprot	Q07869
GeneID	5465;
Target Species	human
Calculated MW	52kd
Concentration	1.0mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000

Note: This product is for in vitro research use only