

Dipeptidyl peptidase 2

Catalog No: #AP82665

Package Size: #AP82665-1 50ug #AP82665-2 100ug #AP82665-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dipeptidyl peptidase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Dipeptidyl aminopeptidase II,Dipeptidyl peptidase 7,Dipeptidyl peptidase II,DPP II,Quiescent cell proline dipeptidase,DPP2,QPP
Accession No.	Uniprot ID:Q9UHL4
Uniprot	Q9UHL4
GeneID	29952;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DPP7

Note: This product is for in vitro research use only