

Actin antibody

Catalog No: #22955

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	Actin antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 146 and 365 of Human ACTA1
Target Name	Actin
Accession No.	Swiss-Prot:P68133Gene ID:58
Uniprot	P68133
GeneID	58;
Concentration	0.1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 42kd

Western blotting: 1:500-1:3000

Images



Sample (30 ug of whole cell lysate)
 A: Raji
 10% SDS PAGE
 Primary antibody diluted at 1: 1000

Background

The product encoded by this gene belongs to the actin family of proteins, which are highly conserved proteins that play a role in cell motility, structure and integrity. Alpha, beta and gamma actin isoforms have been identified, with alpha actins being a major constituent of the contractile apparatus,

while beta and gamma actins are involved in the regulation of cell motility. This actin is an alpha actin that is found in skeletal muscle. Mutations in this gene cause nemaline myopathy type 3, congenital myopathy with excess of thin myofilaments, congenital myopathy with cores, and congenital myopathy with fiber-type disproportion, diseases that lead to muscle fiber defects. [provided by RefSeq]

Note: This product is for in vitro research use only