

DYRK2/4 (Phospho-Tyr382/264) Antibody

Catalog No: #12498



Package Size: #12498-1 50ul #12498-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	DYRK2/4 (Phospho-Tyr382/264) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	DYRK2/4 (Phospho-Tyr382/264) Antibody detects endogenous levels of DYRK2/4 only when phosphorylated at Tyr382/264
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human DYRK2/4 (Phospho-Tyr382/264)
Target Name	DYRK2/4
Modification	Phospho
Other Names	DYRK2
Accession No.	Swiss-Prot#: Q92630/Q9NR20NCBI Gene ID: 8445/8798
Uniprot	Q92630
GeneID	8445;
Target Species	human
Calculated MW	60kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Note: This product is for in vitro research use only