

Hydroxyindole O-methyltransferase

Catalog No: #AP83295



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Hydroxyindole O-methyltransferase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Hydroxyindole O-methyltransferase,HIOMT
Accession No.	Uniprot ID:P10950
Uniprot	P10950
GeneID	281013;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
gene_full_name:ASMT

Note: This product is for in vitro research use only