

MNK2 (Phospho-Thr249) Antibody

Catalog No: #12517



Package Size: #12517-1 50ul #12517-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	MNK2 (Phospho-Thr249) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	MNK2 (Phospho-Thr249) Antibody detects endogenous levels of MNK2 only when phosphorylated at Thr249
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human MNK2 (Phospho-Thr249)
Target Name	MNK2
Modification	Phospho
Other Names	MKNK2, GPRK7, MNK2, MAPK-Interacting Kinase 2
Accession No.	Swiss-Prot#: Q9HBH9NCBI Gene ID: 2872
Uniprot	Q9HBH9
GeneID	2872;
Target Species	human
Calculated MW	52kd
Concentration	1.0mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Images



Western Blot analysis of lysates of mouse liver, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only