

Tissue alpha-L-fucosidase

Catalog No: #AP83367



Package Size: #AP83367-1 50ug #AP83367-2 100ug #AP83367-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Tissue alpha-L-fucosidase |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Bovine |
| Other Names | Alpha-L-fucosidase I,Alpha-L-fucoside fucohydrolase 1,Alpha-L-fucosidase 1 |
| Accession No. | Uniprot ID:Q2KIM0 |
| Uniprot | Q2KIM0 |
| GeneID | 509522; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FUCA1

Note: This product is for in vitro research use only