

Growth hormone secretagogue receptor type 1

Catalog No: #AP83765

Package Size: #AP83765-1 50ug #AP83765-2 100ug #AP83765-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Growth hormone secretagogue receptor type 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	GHS-R,GH-releasing peptide receptor,GHRP,Ghrelin receptor
Accession No.	Uniprot ID:Q95254
Uniprot	Q95254
GeneID	397372;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:GHSR

Note: This product is for in vitro research use only