

## Leptin Receptor (Phospho-Tyr1141) Antibody

Catalog No: #12562

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Leptin Receptor (Phospho-Tyr1141) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu Ms
Specificity	Leptin Receptor (Phospho-Tyr1141) Antibody detects endogenous levels of Leptin Receptor only when phosphorylated at Tyr1141
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human Leptin Receptor (Phospho-Tyr1141)
Target Name	Leptin Receptor
Modification	Phospho
Other Names	LEPR, CD295 antigen, HuB219, LEP-R, OB receptor, LEPRD, OB-R, OBR, CD295, DB, Leptin receptor
Accession No.	Swiss-Prot#: P48357NCBI Gene ID: 3953
Uniprot	P48357
GeneID	3953;
Target Species	human
Calculated MW	96kd
Concentration	1.0mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000

Note: This product is for in vitro research use only