

Major histocompatibility complex class I-related gene protein

Catalog No: #AP83944

Package Size: #AP83944-1 50ug #AP83944-2 100ug #AP83944-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Major histocompatibility complex class I-related gene protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MHC class I-related gene protein,Class I histocompatibility antigen-like protein
Accession No.	Uniprot ID:Q95460
Uniprot	Q95460
GeneID	3140;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MR1

Note: This product is for in vitro research use only