

Serine protease inhibitor Kazal-type 1

Catalog No: #AP84122



Package Size: #AP84122-1 50ug #AP84122-2 100ug #AP84122-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Description

Product Name	Serine protease inhibitor Kazal-type 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Pancreatic secretory trypsin inhibitor,Tumor-associated trypsin inhibitor,TATI,PSTI
Accession No.	Uniprot ID:P00995
Uniprot	P00995
GeneID	6690;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:SPINK1

Note: This product is for in vitro research use only