

Hypoxia-inducible factor 3-alpha

Catalog No: #AP84320

Package Size: #AP84320-1 50ug #AP84320-2 100ug #AP84320-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Hypoxia-inducible factor 3-alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	HIF-3-alpha,Basic-helix-loop-helix-PAS protein MOP7,HIF3-alpha-1,Inhibitory PAS domain protein,IPAS,Member of PAS protein 7,Neonatal and embryonic PAS protein,Mop7,Nepas
Accession No.	Uniprot ID:Q0VBL6
Uniprot	Q0VBL6
GeneID	53417;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hif3a

Note: This product is for in vitro research use only