

## SGK1 (Phospho-Thr256) Antibody

Catalog No: #12667



Package Size: #12667-1 50ul #12667-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	SGK1 (Phospho-Thr256) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu
Specificity	SGK1 (Phospho-Thr256) Antibody detects endogenous levels of SGK1 only when phosphorylated at Thr256
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human SGK1 (Phospho-Thr256)
Target Name	SGK1
Modification	Phospho
Other Names	SGK1, Hsgk, Hsgk1, SGK, Sgk1 variant i3
Accession No.	Swiss-Prot#: O00141NCBI Gene ID: 6446
Uniprot	O00141
GeneID	6446;
Target Species	human
Calculated MW	54kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000

Note: This product is for in vitro research use only