

Low-density lipoprotein receptor-related protein 2

Catalog No: #AP85004

Package Size: #AP85004-1 50ug #AP85004-2 100ug #AP85004-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Low-density lipoprotein receptor-related protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	LRP-2,Glycoprotein 330,gp330,Megalin
Accession No.	Uniprot ID:P98164
Uniprot	P98164
GeneID	4036;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LRP2

Note: This product is for in vitro research use only