

Lupus La protein homolog

Catalog No: #AP85078

Package Size: #AP85078-1 50ug #AP85078-2 100ug #AP85078-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Lupus La protein homolog
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	La autoantigen homolog,La ribonucleoprotein,Ss-b
Accession No.	Uniprot ID:P32067
Uniprot	P32067
GeneID	20823;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ssb

Note: This product is for in vitro research use only