

Ragulator complex protein LAMTOR5

Catalog No: #AP85322



Package Size: #AP85322-1 50ug #AP85322-2 100ug #AP85322-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Description

Product Name	Ragulator complex protein LAMTOR5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	HBV X-interacting protein homolog,HBX-interacting protein homolog
Accession No.	Uniprot ID:Q3SZ68
Uniprot	Q3SZ68
GeneID	516090;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:LAMTOR5

Note: This product is for in vitro research use only