

Metallothionein-1A

Catalog No: #AP85506

Package Size: #AP85506-1 50ug #AP85506-2 100ug #AP85506-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Metallothionein-1A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	MT-1A,Metallothionein-IA,MT-IA,MT1
Accession No.	Uniprot ID:P49068
Uniprot	P49068
GeneID	397417;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MT1A

Note: This product is for in vitro research use only