

PWWP domain-containing protein MUM1

Catalog No: #AP85543



Package Size: #AP85543-1 50ug #AP85543-2 100ug #AP85543-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	PWWP domain-containing protein MUM1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Mutated melanoma-associated antigen 1,MUM-1
Accession No.	Uniprot ID:B1H224
Uniprot	B1H224
GeneID	362838;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Mum1

Note: This product is for in vitro research use only