

Neuronal calcium sensor 1

Catalog No: #AP85641



Package Size: #AP85641-1 50ug #AP85641-2 100ug #AP85641-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Neuronal calcium sensor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	NCS-1,Frequenin homolog,Freq
Accession No.	Uniprot ID:Q8BNY6
Uniprot	Q8BNY6
GeneID	14299;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ncs1

Note: This product is for in vitro research use only