

GFAP (Phospho-Ser13) Antibody

Catalog No: #12757



Package Size: #12757-1 50ul #12757-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	GFAP (Phospho-Ser13) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu
Specificity	GFAP (Phospho-Ser13) Antibody detects endogenous levels of GFAP only when phosphorylated at Ser13
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human GFAP (Phospho-Ser13)
Target Name	GFAP
Modification	Phospho
Other Names	GFAP
Accession No.	Swiss-Prot#: P14136NCBI Gene ID: 2670
Uniprot	P14136
GeneID	2670;
Target Species	human
Calculated MW	50kd
Concentration	1.0mg/ml
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Images



Western Blot analysis of lysates of mouse brain tissue, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only