



Diphosphoinositol polyphosphate phosphohydrolase 3-beta

Catalog No: #AP86272

Package Size: #AP86272-1 50ug #AP86272-2 100ug #AP86272-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Diphosphoinositol polyphosphate phosphohydrolase 3-beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	DIPP-3-beta,Diadenosine 5',5'''-P1,P6-hexaphosphate hydrolase 3-beta,Diadenosine hexaphosphate hydrolase (AMP-forming),Nucleoside diphosphate-linked moiety X motif 11,Nudix motif 11,hAps1,APS1,DIPP3B
Accession No.	Uniprot ID:Q96G61
Uniprot	Q96G61
GeneID	55190;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NUDT11

Note: This product is for in vitro research use only