

CaMKIbeta Antibody

Catalog No: #44301

Package Size: #44301-1 50ul #44301-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CaMKIbeta Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	CaMKIbeta Antibody detects endogenous levels of total CaMKIbeta
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human CaMKIbeta
Target Name	CaMKIbeta
Other Names	PNCK, BSTK3, CaM kinase I beta, CaM kinase IB, CaM kinase Ibeta1, CaM-KI beta, CaMK1b, CaMKI-beta, CaM kinase Ibeta2, CaMKIbeta2
Accession No.	Swiss-Prot#: Q6P2M8NCBI Gene ID:139728
Uniprot	Q6P2M8
GeneID	139728;
Calculated MW	38kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Images



Western blot analysis CaMKIbeta using LOVO whole cell lysates

Note: This product is for in vitro research use only