

Phosphoenolpyruvate carboxykinase, cytosolic [GTP]

Catalog No: #AP86444

Package Size: #AP86444-1 50ug #AP86444-2 100ug #AP86444-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phosphoenolpyruvate carboxykinase, cytosolic [GTP]
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	PEPCK-C
Accession No.	Uniprot ID:P07379
Uniprot	P07379
GeneID	362282;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Pck1

Note: This product is for in vitro research use only