

P2Y purinoceptor 12

Catalog No: #AP86457

Package Size: #AP86457-1 50ug #AP86457-2 100ug #AP86457-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	P2Y purinoceptor 12
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	P2Y12
Accession No.	Uniprot ID:Q9CPV9
Uniprot	Q9CPV9
GeneID	70839;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: P2ry12

Note: This product is for in vitro research use only