

Protein O-mannosyl-transferase 1

Catalog No: #AP86498

Package Size: #AP86498-1 50ug #AP86498-2 100ug #AP86498-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein O-mannosyl-transferase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Dolichyl-phosphate-mannose--protein mannosyltransferase 1
Accession No.	Uniprot ID:Q8R2R1
Uniprot	Q8R2R1
GeneID	99011;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Pomt1

Note: This product is for in vitro research use only