

GBP4 Antibody

Catalog No: #44329



Package Size: #44329-1 50ul #44329-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

| | |
|-----------------------|--|
| Product Name | GBP4 Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | Immunogen affinity purified |
| Applications | WB |
| Species Reactivity | Hu Ms Rt |
| Specificity | GBP4 Antibody detects endogenous levels of total GBP4 |
| Immunogen Type | Peptide |
| Immunogen Description | A synthesized peptide derived from human GBP4 |
| Target Name | GBP4 |
| Other Names | GBP4, GBP-4, Guanylate-binding protein 4, Mpa2, GTP-binding protein 4, Guanylate binding protein 4 |
| Accession No. | Swiss-Prot#: Q96PP9NCBI Gene ID:115361 |
| Uniprot | Q96PP9 |
| GeneID | 115361; |
| Calculated MW | 73kDa |
| Concentration | 1.0mg/mL |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Store at -20°C |


Application Details

Western blotting: 1:1000

Images



Western Blot analysis of RAW cells using GBP4 Polyclonal Antibody



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma, using GBP4 Antibody. The picture on the right is blocked with the synthesized peptide.

Note: This product is for in vitro research use only