

Protein O-linked-mannose beta-1,2-N-acetylglucosaminyltransferase 1

Catalog No: #AP86711

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP86711-1 50ug #AP86711-2 100ug #AP86711-3 1mg

Description

Product Name	Protein O-linked-mannose beta-1,2-N-acetylglucosaminyltransferase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	POMGnT1,UDP-GlcNAc:alpha-D-mannoside beta-1,2-N-acetylglucosaminyltransferase I.2,GnT I.2, MGAT1.2, UNQ746/PRO1475
Accession No.	Uniprot ID:Q8WZA1
Uniprot	Q8WZA1
GeneID	55624;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: POMGNT1

Note: This product is for in vitro research use only