

Ornithine carbamoyltransferase, mitochondrial

Catalog No: #AP86731

Package Size: #AP86731-1 50ug #AP86731-2 100ug #AP86731-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ornithine carbamoyltransferase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Ornithine transcarbamylase,OTCase
Accession No.	Uniprot ID:P11725
Uniprot	P11725
GeneID	18416;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Otc

Note: This product is for in vitro research use only