

IL-1R9 Antibody

Catalog No: #44346

Package Size: #44346-1 50ul #44346-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	IL-1R9 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	IL-1R9 Antibody detects endogenous levels of total IL-1R9
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human IL-1R9
Target Name	IL-1R9
Other Names	IL1RAPL2, IL-1R9, IL-1RAPL-2, IL1RAPL-2, Interleukin 1 receptor 9, IL1RAPL-2-related protein, IL-1-RAPL-2, IL-1R-9, IL1R9, TIGIRR-1, Interleukin-1 receptor 9
Accession No.	Swiss-Prot#: Q9NP60NCBI Gene ID:26280
Uniprot	Q9NP60
GeneID	26280;
Calculated MW	79kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Images



Western blot analysis IL-1R9 using COLO205 whole cell lysates

Note: This product is for in vitro research use only