

Phosphatidylinositol-glycan-specific phospholipase D

Catalog No: #AP87059

Package Size: #AP87059-1 50ug #AP87059-2 100ug #AP87059-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phosphatidylinositol-glycan-specific phospholipase D
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	PI-G PLD, Glycoprotein phospholipase D, Glycosyl-phosphatidylinositol-specific phospholipase D, GPI-PLD
Accession No.	Uniprot ID:O70362
Uniprot	O70362
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Gpld1

Note: This product is for in vitro research use only