

# Olfactory receptor 2C3 Antibody

Catalog No: #44379



Package Size: #44379-1 50ul #44379-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	Olfactory receptor 2C3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	Olfactory receptor 2C3 Antibody detects endogenous levels of total Olfactory receptor 2C3
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human Olfactory receptor 2C3
Target Name	Olfactory receptor 2C3
Other Names	OR2C3, Olfactory receptor 2C4, OR2C4, OR2C5P, OST742, Olfactory receptor OR1-30, Olfactory receptor 2C3, Olfactory receptor 2C5
Accession No.	Swiss-Prot#: Q8N628NCBI Gene ID:81472
Uniprot	Q8N628
GeneID	81472;
Calculated MW	36kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000

## Images



Western blot analysis of lysates from COLO cells, using OR2C3 Antibody. The lane on the right is blocked with the synthesized peptide.

---

Note: This product is for in vitro research use only