

Mothers against decapentaplegic homolog 3

Catalog No: #AP87419



Package Size: #AP87419-1 50ug #AP87419-2 100ug #AP87419-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Description

Product Name	Mothers against decapentaplegic homolog 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MAD homolog 3,JV15-2,SMAD family member 3,SMAD 3,MADH3
Accession No.	Uniprot ID:P84022
Uniprot	P84022
GeneID	4088;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:SMAD3

Note: This product is for in vitro research use only