

ZNF262 Antibody

Catalog No: #44437



Package Size: #44437-1 50ul #44437-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ZNF262 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	ZNF262 Antibody detects endogenous levels of total ZNF262
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human ZNF262
Target Name	ZNF262
Other Names	ZMYM4, CDIR, Cell death inhibiting RNA, MYM, Zinc finger MYM-type protein 4, Zinc finger protein 262, ZNF198L3, Zinc finger, MYM-type 4, ZNF262, KIAA0425
Accession No.	Swiss-Prot#: Q5VZL5NCBI Gene ID:9202
Uniprot	Q5VZL5
GeneID	9202;
Calculated MW	175kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Images



Western Blot analysis of HEK293 lysates, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only