

Protein S100-A11

Catalog No: #AP87987

Package Size: #AP87987-1 50ug #AP87987-2 100ug #AP87987-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein S100-A11
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Calgizzarin, S100 calcium-binding protein A11
Accession No.	Uniprot ID:Q6B345
Uniprot	Q6B345
GeneID	445415;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: S100a11

Note: This product is for in vitro research use only