

TRMT11 Antibody

Catalog No: #44487

Package Size: #44487-1 50ul #44487-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TRMT11 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	TRMT11 Antibody detects endogenous levels of total TRMT11
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human TRMT11
Target Name	TRMT11
Other Names	C6orf75
Accession No.	Swiss-Prot#: Q7Z4G4NCBI Gene ID:60487
Uniprot	Q7Z4G4
GeneID	60487;
Calculated MW	53 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis TRMT11 using 293 whole cell lysates

Note: This product is for in vitro research use only