

# TUSC4 Antibody

Catalog No: #44493



Package Size: #44493-1 50ul #44493-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	TUSC4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	TUSC4 Antibody detects endogenous levels of total TUSC4
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human TUSC4
Target Name	TUSC4
Other Names	NPRL2, G21 protein, Gene 21 protein, NPR2-like protein, NPR2L, TUSC4, Tumor suppressor candidate 4
Accession No.	Swiss-Prot#: Q8WTW4NCBI Gene ID:10641
Uniprot	Q8WTW4
GeneID	10641;
Calculated MW	44 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-3000

## Images



Western Blot analysis of extracts from Jurkat, K562 cells, using NPRL2 Polyclonal Antibody. Secondary antibody was diluted at 1:20000

---

Note: This product is for in vitro research use only