

# Mothers against decapentaplegic homolog 6

Catalog No: #AP88332



Package Size: #AP88332-1 50ug #AP88332-2 100ug #AP88332-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mothers against decapentaplegic homolog 6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	MAD homolog 6,Mad homolog 7,SMAD family member 6,SMAD 6,Madh6,Madh7,Msmad6
Accession No.	Uniprot ID:O35182
Uniprot	O35182
GeneID	17130;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Smad6

Note: This product is for in vitro research use only