

COX82 Antibody

Catalog No: #44511



Package Size: #44511-1 50ul #44511-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	COX82 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	COX82 Antibody detects endogenous levels of total COX82
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human COX82
Target Name	COX82
Other Names	COX8, COX8L
Accession No.	Swiss-Prot#: P10176NCBI Gene ID:1351
Uniprot	P10176
GeneID	1351;
Calculated MW	8 kDa
Concentration	1.0mg/mL
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C

Application Details

WB (1/500 - 1/1000)

Images



Western blot analysis of COX8A expression in A549 (A), SP2/0 (B), rat brain (C) whole cell lysates.

Note: This product is for in vitro research use only