

IL20L2 Antibody

Catalog No: #44516

Package Size: #44516-1 50ul #44516-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	IL20L2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	IL20L2 Antibody detects endogenous levels of total IL20L2
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human IL20L2
Target Name	IL20L2
Other Names	ISG20L2
Accession No.	Swiss-Prot#: Q9H9L3NCBI Gene ID:81875
Uniprot	Q9H9L3
GeneID	81875;
Calculated MW	39 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Note: This product is for in vitro research use only