

# Tumor necrosis factor receptor superfamily member 8

Catalog No: #AP88619



Package Size: #AP88619-1 50ug #AP88619-2 100ug #AP88619-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Tumor necrosis factor receptor superfamily member 8
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	CD30L receptor,Lymphocyte activation antigen CD30,CD30
Accession No.	Uniprot ID:Q60846
Uniprot	Q60846
GeneID	21941;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Tnfrsf8

Note: This product is for in vitro research use only