

Truncated steroid 21-hydroxylase Antibody

Catalog No: #44535

Package Size: #44535-1 50ul #44535-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Truncated steroid 21-hydroxylase Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	Truncated steroid 21-hydroxylase Antibody detects endogenous levels of total Truncated steroid 21-hydroxylase
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human Truncated steroid 21-hydroxylase
Target Name	Truncated steroid 21-hydroxylase
Other Names	Truncated steroid 21-hydroxylase
Accession No.	Swiss-Prot#: Q8TE42NCBI Gene ID:281741
Uniprot	Q8TE42
GeneID	281741;
Calculated MW	8 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis Truncated steroid 21-hydroxylase using 293 whole cell lysates

Note: This product is for in vitro research use only