

Trimethyllysine dioxygenase, mitochondrial

Catalog No: #AP88661

Package Size: #AP88661-1 50ug #AP88661-2 100ug #AP88661-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Trimethyllysine dioxygenase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Epsilon-trimethyllysine 2-oxoglutarate dioxygenase,Epsilon-trimethyllysine hydroxylase,TML hydroxylase,TML-alpha-ketoglutarate dioxygenase,TML dioxygenase,TMLH
Accession No.	Uniprot ID:Q9NVH6
Uniprot	Q9NVH6
GeneID	55217;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TMLHE

Note: This product is for in vitro research use only