

PCTK2 Antibody

Catalog No: #44539

Package Size: #44539-1 50ul #44539-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PCTK2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	PCTK2 Antibody detects endogenous levels of total PCTK2
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human PCTK2
Target Name	PCTK2
Other Names	CDK17, Cyclin-dependent kinase 17, PCTK2, PCTAIRE protein kinase 2, PCTAIRE-motif protein kinase 2, PCTAIRE2, PCTAIRE-2
Accession No.	Swiss-Prot#: Q00537NCBI Gene ID:5128
Uniprot	Q00537
GeneID	5128;
Calculated MW	60 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis of PCTK2 Antibody. The lane on the right is blocked with the PCTK2 peptide.

Note: This product is for in vitro research use only