

TSSK1 Antibody

Catalog No: #44543

Package Size: #44543-1 50ul #44543-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TSSK1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	TSSK1 Antibody detects endogenous levels of total TSSK1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human TSSK1
Target Name	TSSK1
Other Names	SPOGA1, SPOGA4, STK22A, STK22D, TSSK1
Accession No.	Swiss-Prot#: Q9BXA7NCBI Gene ID:83942
Uniprot	Q9BXA7
GeneID	83942;
Calculated MW	42 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Note: This product is for in vitro research use only