

CD300LB Antibody

Catalog No: #44549

Package Size: #44549-1 50ul #44549-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CD300LB Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	CD300LB Antibody detects endogenous levels of total CD300LB
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human CD300LB
Target Name	CD300LB
Other Names	CD300LB, CD300b, CD300b antigen, CMRF35-like molecule 7, CLM7, CLM-7, CMRF35A2, IREM3, IREM-3, LMIR5, TREM5, CMRF35-A2, TREM-5
Accession No.	Swiss-Prot#: A8K4G0NCBI Gene ID:124599
Uniprot	A8K4G0
GeneID	124599;
Calculated MW	23 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis CD300LB using HT-29 whole cell lysates

Note: This product is for in vitro research use only