

E3 ubiquitin-protein ligase TRIM63

Catalog No: #AP88955

Package Size: #AP88955-1 50ug #AP88955-2 100ug #AP88955-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	E3 ubiquitin-protein ligase TRIM63
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Iris RING finger protein,Muscle-specific RING finger protein 1,MuRF-1,RING finger protein 28,RING-type E3 ubiquitin transferase TRIM63,Striated muscle RING zinc finger protein,Tripartite motif-containing protein 63,IRF,MURF1,RNF28,SMRZ
Accession No.	Uniprot ID:Q969Q1
Uniprot	Q969Q1
GeneID	84676;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TRIM63

Note: This product is for in vitro research use only