

CREG1 Antibody

Catalog No: #44569



Package Size: #44569-1 50ul #44569-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	CREG1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	CREG1 Antibody detects endogenous levels of total CREG1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human CREG1
Target Name	CREG1
Other Names	CREG1, CREG, Protein CREG1
Accession No.	Swiss-Prot#: O75629NCBI Gene ID:8804
Uniprot	O75629
GeneID	8804;
Calculated MW	24 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis of lysates from U2OS cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only