

BAI3 Antibody

Catalog No: #44572



Package Size: #44572-1 50ul #44572-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	BAI3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	BAI3 Antibody detects endogenous levels of total BAI3
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human BAI3
Target Name	BAI3
Other Names	BAI3, KIAA0550
Accession No.	Swiss-Prot#: O60242NCBI Gene ID:577
Uniprot	O60242
GeneID	577;
Calculated MW	172 kDa
Concentration	1.0mg/mL
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C

Application Details

WB (1/500 - 1/1000)

Images



Western blot analysis of BAI3 expression in HEK293T (A), Raw264.7 (B), PC12 (C) whole cell lysates.

Note: This product is for in vitro research use only