

V-type proton ATPase catalytic subunit A

Catalog No: #AP89267



Package Size: #AP89267-1 50ug #AP89267-2 100ug #AP89267-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	V-type proton ATPase catalytic subunit A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	V-ATPase 116 kDa isoform a2,Lysosomal H(+)-transporting ATPase V0 subunit a2,TJ6,Vacuolar proton translocating ATPase 116 kDa subunit a isoform 2
Accession No.	Uniprot ID:Q9Y487
Uniprot	Q9Y487
GeneID	23545;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ATP6V1A

Note: This product is for in vitro research use only