

Enoyl-CoA hydratase, mitochondrial

Catalog No: #AP89383

Package Size: #AP89383-1 50ug #AP89383-2 100ug #AP89383-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Enoyl-CoA hydratase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Enoyl-CoA hydratase 1,Short-chain enoyl-CoA hydratase,SCEH
Accession No.	Uniprot ID:P30084
Uniprot	P30084
GeneID	1892;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ECHS1

Note: This product is for in vitro research use only